

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	32606-5K	2026-10-31
Low	32606-51	2026-10-31
Normal	32606-52	2026-11-04
High	32606-53	2026-11-06

MEDONIC

		CONTROL	L
		LOT	32606-51
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	67.8	7.0
RDW	%	15.6	4.5
Hct	%	15.1	3.0
Plt	10 ⁹ /L	82	15
MPV	fL	11.2	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.7	2.0
MCHC	g/dL	36.4	4.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	44.1	7.0
Mid	%	5.2	5.0
Gran	%	50.7	7.0

		CONTROL	N
		LOT	32606-52
		Mean ²	± Range ³
		4.18	0.18
		83.6	7.0
		12.9	4.5
		34.9	3.5
		235	30
		11.1	1.5
		8.6	0.6
		12.3	0.4
		29.4	2.0
		35.2	4.0
		4.0	0.6
		0.5	0.5
		4.1	1.0
		46.4	5.0
		6.1	5.0
		47.5	5.0

		CONTROL	H
		LOT	32606-53
		Mean ²	± Range ³
		4.98	0.22
		92.1	7.0
		13.2	4.5
		45.9	4.0
		529	60
		11.5	1.5
		20.4	1.8
		15.9	0.5
		31.9	2.0
		34.6	4.0
		12.5	1.5
		1.2	1.0
		6.7	2.5
		61.3	5.0
		5.7	5.0
		33.0	5.0

1
bg Трайност на отворен флакон 14 дни (bg)
hr Stabilität otvorene bočice 14 dni (hr)
cs Stabilita otevřené ampulky 14 dní (cs)
da Åben rør stabilitet 14 dage (da)
et Stabiilsus avatud vialis 14 päevade (et)
fr Stabilité en flacon ouvert 14 jour (fr)
de Stabilität geöffneter Flaschen 14 tage (de)
el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it Stabilità della fiala aperta 14 giorni (it)
lv Atvērt flakonu stabilitāti 14 dienas (lv)
lt Stabilitumas atidarius buteliuką 14 dienas (lt)
no Åpen rør/glass stabilitet 14 dager (no)
pl Trwałość po otwarciu fiolek 14 dni (pl)
pt Estabilidade após abertura do frasco 14 dias (pt)
ro Termenul de valabilitate al fiolei deschise 14 zile (ro)
ru Стабильность вскрытого флакона 14 дней (ru)
sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es Estabilidad de la cápsula abierta 14 días (es)
sv Hållbarhet för öppnad flaska 14 dagar (sv)
tr Açık şişe stabilitesi 14 gün (tr)

2
Sredno
Sredstvo
Průměr
Middelverdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3
Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı



www.boule.com

Boule



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25122-5

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 32606-51
Prod.: Con-Diff Low
ID: 3260651+
Exp.: 31-Oct-2026

Instr.: M-series
Lot: 32606-52
Prod.: Con-Diff Norm
ID: 3260652+
Exp.: 4-Nov-2026

Instr.: M-series
Lot: 32606-53
Prod.: Con-Diff High
ID: 3260653+
Exp.: 6-Nov-2026

